



Next-Generation Protein Sequencing on Quantum-Si's Platinum™: Discernment of Protein Mixtures and Peptide Modifications


Kenneth A. Skinner, Haidong Huang, Michael J. Warchol, Mathivanan Chinnaraj, Khanh D. Q. Nguyen, Ben Moree, Katherine Johnson, and Brian Reed

Quantum-Si Incorporated, 29 Business Park Dr., Branford, CT 06405

INTRODUCTION


Peptide-level analysis coupled with amino acid resolution is important to elucidate how proteins influence health and disease. Peptide-centric methodologies using liquid chromatography-tandem mass spectrometry (LC-MS/MS) are widely used to detect native and modified proteins. While LC-MS/MS stands as the dominant bottom-up proteomics technique, accessing LC-MS/MS and interpreting MS/MS spectra remains challenging. Due to the large capital expenditure and space required for LC-MS/MS, core facilities are often required for processing and analyzing protein samples. In addition, unexpected fragmentation patterns can preclude peptide detection, requiring additional interrogation of peptides for a more complete analysis. To address these hurdles, we demonstrate the ability of Quantum-Si's Platinum™ next-generation benchtop protein sequencer to identify proteins in mixtures and discern endogenous modifications from instrument-related effects with single molecule resolution, complementing existing LC-MS/MS techniques. This dynamic protein sequencing approach employs a mixture of dye-labeled N-terminal amino acid recognizers (NAARs) and aminopeptidases to probe digested peptides. The order of recognizer binding and kinetic properties of recognition segments are analyzed to determine peptide sequence and associated proteins.

Here, we identify protein mixtures consisting of therapeutically relevant growth factors, cytokines, and secreted proteins. Moreover, kinetic signatures from NAARs not only reveal peptides that escape MS/MS mapping due to factors such as peptide length and pyroglutamate formation, but also enable differentiation between Asparagine (Asn) and deamidated peptide variants. These findings underscore Platinum as an alternative and complementary technique to LC-MS/MS for protein identification, peptide-level mapping, and monitoring critical quality attributes (CQAs) during product development.



In proteomics analysis, the concern is always *what is missed*

METHODS



Platinum


RESULTS AND DISCUSSION

Platinum reveals protein components in mixtures




Green boxes indicate peptides sequenced on Platinum.

Platinum detects a short FGF2 peptide that escaped MS matching




Deamidation alters the structure and function of biomolecules




REFERENCE

Brian D. Reed et al. Science, 2022, 378 (6166) 186–192.
Preeti Purwaha et al. Analytical Chemistry, 2014, 86 (12) 5633-5637.
Josef Vlasek et al. Analytical Biochemistry, 2009, 392 (2) 145-154.
Lukas Käll et al. Journal of Proteome Research, 2008, 7 (01) 40–44.
Kai-Fa Huang et al. PNAS, 2005, 102 (37) 13117-13122.

Platinum distinguishes Asn-containing and deamidated CDFN peptides



Platinum sequences an IL6 peptide that escaped MS matching



(A) Gln cyclization forms neutral lactam group. (B) TIC chromatogram of IL6 subjected to bottom-up sample prep. (C) List of calc'd QIRYILDGIALRK fragment ions. (D) Biomolecule fragmentation spectra of recombinant human IL6. (E) Relative TIC ratio of the b3 fragment peaks from synthetic QIRYILDGIALRK depends on Orbitrap HCD collision energy. (F) IL6 sequencing trace.

TRADEMARKS/LICENSING

All trademarks are the property of Quantum-Si, Inc. or their respective owners. For specific trademark information, see www.quantum-si.com/privacy-policy.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

SCAN TO VISIT



QUANTUM-SI.COM